Otto Bock®


QUALITY FOR LIFE

28U22

(D)	Otto Bock Walk On® Flex Unterschenkelorthese	3
(GB)	Otto Bock Walk On® Flex AFO orthosis	6
(F)	Orthèse tibiale Otto Bock Walk On® Flex	9
(I)	Ortesi di gamba Walk On® Flex Otto Bock	12
(E)	Otto Bock Walk On® Flex Órtesis tibial	15
(P)	Ortótese de joelho, tornozelo e pé Walk On® Flex Otto Bock	18
(NL)	Otto Bock Walk On® Flex onderbeenorthese	21
(S)	Otto Bock Walk On® Flex underbensortos	24
(DK)	Otto Bock Walk On® Flex underekstremitetsortose	27
(N)	Otto Bock Walk On® Flex skinnleggsortese	30
(FIN)	Otto Bock Walk On® Flex-sääriortoosi	33
(PL)	Orteza podudzia Walk On® Flex firmy Otto Bock	36
(CZ)	Kotníková ortéza (AFO) Otto Bock Walk On® Flex	39
(TR)	Otto Bock Walk On® Flex alt ekstremite ortezi	42
(GR)	Ποδοκνημικός κηδεμόνας Walk On® Flex της Otto Bock	45
(RUS)	Ортез голени Walk On® Flex производства компании Otto Bock	48
(J)	ウォーク オン フレックス 短下肢装具	51

	Artikelnr./Article No.		Schuhgröße/ Shoe Size
	links/left	rechts/right	
	28U22=L small	28U22=R small	36 - 39
	28U23=L medium	28U23=R medium	39 - 42
	28U23=L large	28U23=R large	42 - 45
	28U23=L X-large	28U23=R X-large	45 - 48

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⚠ WARNING	Warnings regarding possible risks of severe accident or injury.
⚠ CAUTION	Warnings regarding possible risks of accident or injury.
NOTICE	Warnings regarding possible technical damage.
INFORMATION	Additional information on the fitting / use.

INFORMATION

Please read the following directions carefully and pay special attention to the safety information!

1 Description

1.1 Intended Use

28U22 Walk On® Flex is to be used **exclusively** for orthotic treatment of the lower extremities and is exclusively suitable for contact with healthy skin.

1.2 Indications and Effects

Mild to moderate drop foot (weakness of dorsiflexion) with no or only mild spasticity, e.g. after peroneal palsy or stroke.

The orthosis is designed to help people with above described problems to achieve a more natural gait pattern. The Walk On® Flex lifts the foot during swing phase and limits plantar flexion.

NOTICE

The Walk On® Flex is a dynamic ankle foot orthosis and must not be used for prevention of contractures etc.

Indication must be determined by the physician.

1.3 Contraindication

1.3.1 Relative contraindication

In case of the following indications, consultation with a physician is required:

Skin diseases/-injuries, inflammatory symptoms, prominent scars that are swollen, reddened and feel excessively warm when touched; Lymphatic flow disorders, including unclear soft tissue swelling distant to the body area to which the appliance will be applied; Sensory and circulatory disorders of the lower extremity.

1.4 Safety Instructions

INFORMATION

The patient is to be instructed in the proper use/care of the product.

The initial fitting and application of the product must be carried out by trained, qualified personnel only. The daily duration of use and period of application are dependent on medical indication by the physician.

⚠ WARNING

Risk of accident when driving a motor vehicle. The ability of orthosis users to drive a vehicle is determined on a case-by-case basis. Criteria include the type of fitting (clinical picture, fitting, design of the orthosis) and the individual abilities of the orthosis user.

All persons are required to observe their country's national and state driving laws when operating motor vehicles. For insurance purposes, drivers should have their driving ability examined and approved by an authorised test centre.

⚠ CAUTION

Risk of injury as a result of improper use. The product is designed for use on **one** patient only. Consult a physician immediately if you experience unusual changes (such as increase in pain).

Improper changes to the product are not permitted.

An orthosis/support applied too tightly to the body can cause local pressure and, in some cases, even restrict adjacent vessels or nerves. Do not apply the product too tightly.

⚠ CAUTION

Risk of injury in exceptional situations. The user **must be informed of the risks** that exceptional situations might present. For example, jumping down from great height may expose the spring to severe overload and cause it to break; or the orthosis could prevent the wearer from braking quickly and suddenly while driving a car.

NOTICE

Damage due to incorrect environmental conditions. This product is not flame-resistant. Keep the product away from flames or other heat sources.

The product should not come into contact with grease or acidic agents, unguents and lotions.

1.5 Construction

Thanks to its construction with three dynamic elements (the pylon, the heel and forefoot of the sole) the orthosis provides a smooth natural gait even on uneven ground.

NOTICE

Damage due to inadmissible handling. The orthosis is made of pre-impregnated glass fiber material and cannot be thermoformed. The foot sole and the connection elements must be free from holes as these would interrupt the fibers and weaken the component.

2 Handling

2.1 Size Selection

The size selection of the orthosis is based on the shoe size. Left and right. S–XL.

2.2 Adapting and Applying the product

For adjustments, the sole can be cut with a pair of scissors.

Adjust the sole:

- 1) If the patient has a removable sole in his or hers shoes, use that for marking the correct size at the Walk On Flex trimable sole (Fig. 1).
- 2) Cut the trimable sole after your markings (Fig. 2), the lateral side should not be cut more than is required, so that the medial part of the spring does not press against the medial side of the foot including the medial malleolus (Fig. 3). If the sole requires more reduction in width, then cut bilaterally.

2.3 Instruction for use and care

Material: Glass fiber composite, PCM®, Velcro®.

Cleaning: Soft calf band: Machine wash, 40° with a standard mild detergent. Wash out thoroughly and air-dry. Note: Residues of cleaning agents may cause skin irritations and wear of material.

Upper limit temperature: 120° C.

Disposal after use: Combustible material.

3 Further Usage Restrictions

The product has been designed for use on only **one** patient. The daily duration of use and period of application are dependent on the medical indication.

CAUTION

Danger of injury due to improper use. Parts that are shaped for a particular patient and those parts that come into contact with the skin may present functional or hygienic risks when used by another person.

Latex Safe: To our knowledge this product does not contain any natural rubber.

4 Liability

The manufacturer's warranty applies only if the device has been used under the conditions and for the purposes described. The manufacturer recommends that the device be used and maintained according to the instructions for use.

5 CE Conformity

This device meets the requirements of the 93/42/EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Otto Bock with sole responsibility according to appendix VII of the guidelines.

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Otto Bock HealthCare GmbH