Otto Bock®

QUALITY FOR LIFE


8356 Genu Direxa

(D)	Knieorthese.....	5
(GB)	Knee Support	8
(F)	Orthèse de genou	11
(I)	Ginocchiera	14
(E)	Rodillera	18
(NL)	Kniebandage	21
(S)	Knäskydd	24
(DK)	Knæortose	27
(FIN)	Polviortoosi	30
(PL)	Orteza kolanowa	33
(CZ)	Kolenní ortéza	37
(HR)	koljena ortoza	43
(SK)	Ortéza kolena	46
(TR)	Diz ortezi	49
(RUS)	Коленный ортез	52
(中国)	膝关节矫形器	56

	Größe/ Size	Umfang/Circumference*			
		D (cm)	D (inch)	F (cm)	F (inch)
	XXS	29–32	11,4–12,6	36–40	14,2–15,8
	XS	32–35	12,6–13,8	40–44	15,8–17,3
	S	35–38	13,8–15,0	44–48	17,3–18,9
	M	38–41	15,0–16,1	48–52	18,9–20,5
	L	41–44	16,1–17,3	52–56	20,5–22,0
	XL	44–48	17,3–18,9	56–61	22,0–24,0
	XXL	48–51	18,9–20,1	61–67	24,0–26,4

* unterhalb/oberhalb Patellamitte / below/above mid patella

Material	PA, EA, PU foam, aluminium
----------	----------------------------

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Explanation of Symbols

English

CAUTION Warnings regarding possible risks of accident or injury.

NOTICE Warnings regarding possible technical damage.

INFORMATION Additional information on the fitting / use.

INFORMATION

Please read these instructions for use carefully and pay special attention to the safety information!

1 Function

Stabilising and relieving knee support. Double-jointed aluminium bars along the sides. User-friendly brace which is anatomically designed. A material which drains away moisture and keeps the skin dry.

2 Indications

Knee pain with light to moderate ligament laxities and/or perception of instability:

- Osteoarthritis with light to moderate degenerative ligament laxities and/or perception of instability (cruciate ligament accentuated).
- Post-traumatic and post-operative ligament laxities or perception of instability.
- Rheumatoid Arthritis.
- Patello-femoral pain syndrome (e.g. chondromalacia, chondropathy of the patella, beginning osteoarthritis of the patellofemoral joint).

2.1 Contraindications

2.1.1 Absolute Contraindications

None.

2.1.2 Relative Contraindications

In case of the following indications, consultation with a physician is required: Skin diseases/injuries, inflammatory symptoms, prominent scars on body parts to which the support is to be applied that are swollen, reddened and feel excessively warm; lymphatic flow disorders, including unclear soft tissue swelling distant to the body area to which the support is to be applied; sensory and circulatory disorders.

3 Safety Instructions

INFORMATION

The patient is to be instructed in the proper use/care of the product.

The initial fitting and application of the product must be carried out by trained, qualified personnel only. The daily duration of use and period of application are dependent on medical indication by the physician.

CAUTION

Risk of injury as a result of improper use. The product is designed for use on one patient only. For hygienic reasons, stockings are always to be worn under the orthosis.

An orthosis/support applied too tightly to the body can cause local pressure and, in some cases, even restrict adjacent vessels or nerves. Do not apply the product too tightly. Consult a physician immediately if you experience unusual changes (such as increase in pain).

Improper changes to the product are not permitted.

CAUTION

Risks due to high temperatures. When an orthosis gets into contact with great heat for a longer time (e.g. longer exposure to cigarette embers), the orthosis' material can melt and cause burns on the patient's skin. Keep away from all sources of ignition, especially lit cigarettes.

NOTICE

Risk of damage due to incorrect environmental conditions. This product is not flame-resistant. Keep the product away from flames or other heat sources.

The product should not come into contact with grease or acidic agents, unguents and lotions, as this can reduce the useful life – especially of pads.

4 How to apply

The bars are set individually by trained personnel.

- 1) Pull on the knee support so that the hole at the front is aligned with the kneecap.
- 2) First of all close the lower and then the upper Velcro® fastener.

5 Washing instructions

- 1) Remove detachable bars.
- 2) Close all Velcro® fasteners.
- 3) Wash at 40 °C / 104 °F using a laundry bag, if possible.
- 4) Rinse thoroughly so that all residual detergent is removed.
- 5) Stretch while wet, and then dry.
- 6) Insert the bars – product number on the bar must be visible from the outside (Fig. 3).

6 Further Usage Restriction

The product is designed for use on one patient only. Re-use of the product is not permitted.

The daily duration of use and period of application are dependent on the medical indication.

7 Liability

The manufacturer's warranty applies only if the device has been used under the conditions and for the purposes described. The manufacturer recommends that the device be used and maintained according to the instructions for use.

8 CE Conformity

This device meets the requirements of the 93/42/EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Otto Bock with sole responsibility according to appendix VII of the guidelines.

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Otto Bock HealthCare GmbH