Otto Bock®

QUALITY FOR LIFE


8357 Genu Direxa Stable

(D)	Knieorthese.....	5
(GB)	Knee Support.....	9
(F)	Orthèse de genou	12
(I)	Ginocchiera	16
(E)	Rodillera.....	20
(NL)	Kniebandage	24
(S)	Knäskydd	28
(DK)	Knæortose	32
(FIN)	Polviortoosi.....	35
(PL)	Orteza kolanowa	39
(SK)	Ortéza kolena.....	43
(TR)	Diz ortezi.....	47
(RUS)	Коленный ортез.....	50
(中国)	膝关节矫形器	55

①

	Größe/ Size	Umfang/Circumference*			
		D (cm)	D (inch)	F (cm)	F (inch)
	XXS	29–32	11,4–12,6	36–40	14,2–15,8
	XS	32–35	12,6–13,8	40–44	15,8–17,3
	S	35–38	13,8–15,0	44–48	17,3–18,9
	M	38–41	15,0–16,1	48–52	18,9–20,5
	L	41–44	16,1–17,3	52–56	20,5–22,0
	XL	44–48	17,3–18,9	56–61	22,0–24,0
	XXL	48–51	18,9–20,1	61–67	24,0–26,4

* unterhalb/oberhalb Patellamitte / below/above mid patella

Material	PA, EA, PU foam, aluminium
----------	----------------------------

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⚠ CAUTION Warnings regarding possible risks of accident or injury.

NOTICE Warnings regarding possible technical damage.

INFORMATION Additional information on the fitting / use.

INFORMATION

Please read these instructions for use carefully and pay special attention to the safety information!

1 Function

Extra-stabilising and relieving knee support. Double-jointed and lockable aluminium bars along the sides. User-friendly brace which is anatomically designed. A material which drains away moisture and keeps the skin dry.

2 Indications

Knee pain with light to moderate ligament laxities and/or perception of instability with indication for permanent or temporary limitation of flexion and/or extension:

- osteoarthritis with light to moderate degenerative ligament laxities and/or perception of instability (ACL, PCL)
- post-traumatic and post-operative ligament laxities or perception of instability
- Rheumatoid Arthritis
- Patello-femoral pain syndrome (e.g. chondromalacia, chondropathy of the patella, beginning osteoarthritis of the patellofemoral joint)

2.1 Contraindications

2.1.1 Absolute Contraindications

None.

2.1.2 Relative Contraindications

In case of the following indications, consultation with a physician is required: Skin diseases/injuries, inflammatory symptoms, prominent scars on body parts to which the support is to be applied that are swollen, reddened and

feel excessively warm; lymphatic flow disorders, including unclear soft tissue swelling distant to the body area to which the support is to be applied; sensory and circulatory disorders.

3 Safety Instructions

INFORMATION

The patient is to be instructed in the proper use/care of the product.

The initial fitting and application of the product must be carried out by trained, qualified personnel only. The daily duration of use and period of application are dependent on medical indication by the physician.

CAUTION

Risk of injury as a result of improper use. The product is designed for use on one patient only. For hygienic reasons, stockings are always to be worn under the orthosis.

An orthosis/support applied too tightly to the body can cause local pressure and, in some cases, even restrict adjacent vessels or nerves. Do not apply the product too tightly. Consult a physician immediately if you experience unusual changes (such as increase in pain).

Improper changes to the product are not permitted.

CAUTION

Risks due to high temperatures. When an orthosis gets into contact with great heat for a longer time (e.g. longer exposure to cigarette embers), the orthosis' material can melt and cause burns on the patient's skin. Keep away from all sources of ignition, especially lit cigarettes.

NOTICE

Risk of damage due to incorrect environmental conditions. This product is not flame-resistant. Keep the product away from flames or other heat sources.

The product should not come into contact with grease or acidic agents, unguents and lotions, as this can reduce the useful life – especially of pads.

4 How to apply

The bars are set individually by trained personnel.

1) Adjusting the Joints using the Click 2 Go System (Fig. 2)

In order to exchange the stops; Insert the desired extension and flexion stops. When inserting the stops they will engage with slight resistance.

Available (functional angle) extension stops:

0°, 10°, 20°, 30°, 45°.

Available (functional angle) flexion stops:

0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°.

Close the Velcro® tape around the joint.

INFORMATION

Both medial and lateral joints must be fitted with the same degree stops!

- 2) Pull up the knee support so that the hole in front is aligned with the kneecap.
- 3) Tighten the band at the back on the calf and thigh. The back band may well be left closed after initial adjustment, even when washing.
- 4) Tighten the band at the front on the calf and thigh.

5 Washing instructions

- 1) Remove detachable bars.
- 2) Close all Velcro® fasteners.
- 3) Wash at 40 °C / 104 °F using a laundry bag, if possible.
- 4) Rinse thoroughly so that all residual detergent is removed.
- 5) Stretch while wet, and then dry.
- 6) Insert the bars – product number on the bar must be visible from the outside (Fig. 3)

6 Further Usage Restriction

The product is designed for use on one patient only. Re-use of the product is not permitted.

The daily duration of use and period of application are dependent on the medical indication.

7 Liability

The manufacturer's warranty applies only if the device has been used under the conditions and for the purposes described. The manufacturer recommends that the device be used and maintained according to the instructions for use.

8 CE Conformity

This device meets the requirements of the 93/42/EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Otto Bock with sole responsibility according to appendix VII of the guidelines.

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Otto Bock HealthCare GmbH