Otto Bock®


QUALITY FOR LIFE

8360 Genu Carezza Patella Stabilizer

(D)	Kniebandage	4
(GB)	Knee Support	7
(F)	Genouillère	10
(I)	Ginocchiera	14
(E)	Rodillera	17
(NL)	Kniebandage	21
(S)	Knäskydd	24
(DK)	Knæbandage	28
(FIN)	Polvituki	31
(PL)	Orteza stawu kolanowego	34
(CZ)	Kolenní bandáž	37
(HR)	zavoj za koljeno	41
(TR)	Dizlik	44
(RUS)	Коленный ортез	47
(中国)	护膝	51


	Größe/ Size	Umfang/Circumference*			
		D (cm)	D (inch)	F (cm)	F (inch)
	XXS	29–32	11,4–12,6	36–40	14,2–15,8
	XS	32–35	12,6–13,8	40–44	15,8–17,3
	S	35–38	13,8–15,0	44–48	17,3–18,9
	M	38–41	15,0–16,1	48–52	18,9–20,5
	L	41–44	16,1–17,3	52–56	20,5–22,0
	XL	44–48	17,3–18,9	56–61	22,0–24,0
	XXL	48–51	18,9–20,1	61–67	24,0–26,4

* unterhalb/oberhalb Patellamitte / below/above mid patella

Material	PA, EA, PU foam
-----------------	-----------------

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Explanation of Symbols

English

⚠ CAUTION Warnings regarding possible risks of accident or injury.

NOTICE Warnings regarding possible technical damage.

INFORMATION Additional information on the fitting / use.

INFORMATION

Please read these instructions for use carefully and pay special attention to the safety information!

1 Function

Stabilises the patella laterally and has a pain-relieving effect. This is a small support with a important functions. Easy to use when moving around and taking part in sports. Can be used on both the right and the left knee. User-friendly support which is anatomically designed. A material which drains away moisture and keeps the skin dry.

2 Indications

- Patellofemoral pain syndrome (e.g. chondromalacia, chondropathia patellae, moderate arthrosis of the patellofemoral joint)
- Patellofemoral maltracking
- Stabilization after patellar luxations/subluxations
- Stabilization after lateral release operation

2.1 Contraindications

2.1.1 Absolute Contraindications

None.

2.1.2 Relative Contraindications

In case of the following indications, consultation with a physician is required: Skin diseases/injuries, inflammatory symptoms, prominent scars on body parts to which the support is to be applied that are swollen, reddened and feel excessively warm; lymphatic flow disorders, including unclear soft tissue swelling distant to the body area to which the support is to be applied; sensory and circulatory disorders.

3 Safety Instructions

INFORMATION

The patient is to be instructed in the proper use/care of the product.

The initial fitting and application of the product must be carried out by trained, qualified personnel only. The daily duration of use and period of application are dependent on medical indication by the physician.

CAUTION

Risk of injury as a result of improper use. The product is designed for use on one patient only. For hygienic reasons, stockings are always to be worn under the orthosis.

An orthosis/support applied too tightly to the body can cause local pressure and, in some cases, even restrict adjacent vessels or nerves. Do not apply the product too tightly. Consult a physician immediately if you experience unusual changes (such as increase in pain).

Improper changes to the product are not permitted.

CAUTION

Risks due to high temperatures. When an orthosis gets into contact with great heat for a longer time (e.g. longer exposure to cigarette embers), the orthosis' material can melt and cause burns on the patient's skin. Keep away from all sources of ignition, especially lit cigarettes.

NOTICE

Risk of damage due to incorrect environmental conditions. This product is not flame-resistant. Keep the product away from flames or other heat sources.

The product should not come into contact with grease or acidic agents, unguents and lotions, as this can reduce the useful life – especially of pads.

4 How to apply

- 1) Pull the brace on so that the opening fits around the patella, the central seam is right at the front and the cushion ends up on the outside of the patella. (In certain cases the cushion is to be placed on the inside, contact your care provider if you feel unsure.)
- 2) Thread the lower band into the upper plastic ring, pull back and do up the Velcro® fastener.
- 3) Thread the upper band into the lower plastic ring, pull back and fasten the Velcro®.
- 4) The band regulates the pressure of the cushion against the patella.
- 5) Now the cushion will lie permanently against the side of the patella.

5 Washing instructions

- 1) Close all Velcro® fasteners.
- 2) Wash at 40 °C / 104 °F using a laundry bag, if possible.
- 3) Rinse thoroughly so that all residual detergent is removed.
- 4) Stretch while wet, and then dry.

6 Further Usage Restriction

The product is designed for use on one patient only. Re-use of the product is not permitted.

The daily duration of use and period of application are dependent on the medical indication.

7 Liability

The manufacturer's warranty applies only if the device has been used under the conditions and for the purposes described. The manufacturer recommends that the device be used and maintained according to the instructions for use.

8 CE Conformity

This device meets the requirements of the 93/42/EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Otto Bock with sole responsibility according to appendix VII of the guidelines.

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Otto Bock HealthCare GmbH