


Otto Bock®

QUALITY FOR LIFE


50R20

(D)	Dorso Osteo Care Rumpforthese	5
(GB)	Dorso Osteo Care Spinal Orthosis	8
(F)	Orthèse de tronc Dorso Osteo Care	12
(I)	Ortesi per busto Dorso Osteo Care	15
(E)	Órtesis de tronco Dorso Osteo Care	19
(P)	Ortótise de tronco Dorso Osteo Care	22
(NL)	Dorso Osteo Care romporthese	26
(S)	Dorso Osteo Care kroppsRTOS	29
(DK)	Dorso Osteo Care kroppsRTose	33
(PL)	Orteza tułowia Dorso Osteo Care	36
(H)	Dorso Osteo Care törzsRTézis	39
(CZ)	Trupová ortéza Dorso Osteo Care	42
(HR)	Ortoza za kralježnicu kod osteoporoz	45
(TR)	Dorso Osteo Care Osteoporoz Ortezi	48
(GR)	Όρθωση κορμού Dorso Osteo Care	52
(RUS)	Ортез туловища Dorso Osteo Care	55
(J)	Dorso Osteo Care 身体矯正具	59
(中國)	Dorso Osteo Care 軀干矯正器	62

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1 Application

English

The 50R20 Dorso Osteo Care Spinal Orthosis is to be **exclusively** used for orthotic fitting of the trunk and is designed for the contact with the intact skin only.

2 Indication and Effects

The 50R20 Dorso Osteo Care Orthosis is indicated for osteoporotic changes in the lumbar and thoracic areas of the spine and follow-up treatment of the osteoporotic vertebral body fractures. The Dorso Osteo Care orthosis can be employed during the stages of osteoporosis in which muscle action can be used and increased. The indication must be diagnosed by the physician.

The upper body is raised into an upright position by retracting the shoulder girdle and abdomen. This effect is achieved by setting straps and adjustable back plate. Thoracic and abdominal breathing are not restricted and the back muscles are activated.

3 Size Selection

The orthosis is available in two sizes and the size selection is based on the trunk length.

Size	for trunk length in cm (approx. 3 cm below C7 to approx. 2 cm above the gluteal fold)
50R20=S	39 – 41.5
50R20=M	42 – 44.5

4 Scope of Delivery and Design

4.1 Scope of Delivery (Fig. 1–3)

- | | |
|-----------------------|---|
| (A) Back Plate | (G) Plastic Buckle |
| (B) Back Pad | (H) Velcro Straps |
| (C) Shoulder Strap | (I) Plastic Loop |
| (D) Waist Strap | (J) Clip Rivets |
| (E) Abdominal Support | (K) Trunk Buckles |
| (F) Velcro Pads | (L) Brace Pad with Double-sided Adhesive Tape |

4.2 Design (Fig. 1–3)

The 50R20 Dorso Osteo Care is supplied ready for trial fitting. The Back Plate (A) consists of coated aluminum and can be easily adjusted. The Back Pad (B) with hook & loop attachment is made of the SpaceTex® spacer fabric, which is breathable, gentle to the skin and washable. The Shoulder Strap (C) is connected to the Back Plate (A) with plastic loops (1–4) and consists of a fleece strap with inserts padded with semi-elastic SpaceTex® spacer fabric in the axillary area. The Shoulder Strap (C) is provided with a safe and simple Plastic Buckle (G), which is located at the front. The Abdominal Support (E) is made of SpaceTex® spacer fabric and is connected to the Back Plate (A) via the Waist Strap (D) with plastic loops (5–8).

5 Adjusting the Orthosis

The initial fitting and application of the 50R20 Dorso Osteo Care spinal orthosis must be carried out by trained and qualified personnel. The Back Plate (A) can be cold-shaped by hand without use of any tools (Fig. 7). To do this, detach the velcroed Back Pad (B), Waist Strap (D) and Shoulder Strap (C) from the Back Plate (A). Adapt the Back Plate (A) exactly to the course of spinal column while observing the specified measuring points (C7 / gluteal fold). Cavities should be avoided, as should any pad pressure. Subsequently, velcro the Back Pad (B) and attach the Shoulder Strap (C) following the numbered sequence (1–4), the Waist Strap (D) with plastic loops (5–8).

6 Applying and Trial Fitting

The orthosis is pulled on following the rucksack principle (Fig. 4). First, open the Plastic Buckle (G) of the Shoulder Strap (C) and the hook & loop closure on one side of the Abdominal Support (right or left, whichever you choose). Pull on the orthosis like a rucksack and position the Back Plate (A) with the Back Pad (B) so that it fits the back. Position the Abdominal Support (E) in the center of the abdomen and close the hook & loop closure (Fig. 5). Then close the Plastic Buckle (G) below the breast (Fig. 6). Adjust the Shoulder Strap (C) and Waist Strap (D) to the proper length using the Velcro Pads. The Shoulder Strap (C) can be additionally trimmed in length. To do so, remove the Plastic Buckle (G) (Fig. 8), cut the Shoulder Strap uniformly on both sides and subsequently re-attach the Plastic Buckle with hook & loop closures. The Plastic Buckle (G) must always be located in the center of the body.

To ensure easier application of the orthosis in patients with limited motility, the enclosed Trunk Buckles (K) can be mounted. To do so, loosen the plastic loops (3+4) of the Shoulder Strap (C) from the Back Plate.

Screw the Trunk Buckles (K) to the Back Plate (A) following the numbered sequence (3+4) (Fig. 9) and glue the enclosed Pads (L) to the Trunk Buckles and cold-shape them by hand to optimally fit the trunk. Hook the plastic loops (3+4) of the Shoulder Strap (C) to the Trunk Buckles (Fig. 10).

To prevent constriction, do not make the Shoulder Strap (C) too tight. It should be possible to put a finger easily between the strap and the shoulder.

The Abdominal Pad (E) positioned in the center must tightly fit the trunk. Set the Waist Strap (D) to the proper length to provide a tight fit. The Shoulder Strap should not be twisted.

7 After the Trial Fitting

Once the final fitting has been completed, replace the Velcro Pads with Clip Rivets (J) (Fig. 11). To do so, rivet the Shoulder and Waist Straps together by hand. Use the large clip rivet for the Shoulder Strap (numbers 1–4) and the Waist Strip (numbers 7+8) and the small rivet for the Waist Strip (numbers 5+6).

Explain to the patient how to apply the orthosis and ensure its correct fit. Practice the donning and doffing of the orthosis with the patient.

8 Side Effects

Attention!

All aids – supports and orthoses – which are applied on the trunk may result in local pressure spots or, in some cases, even restrict the adjacent vessels or nerves when they are applied too tightly.

9 Contraindications

Attention!

In the case of the following symptoms the application and wearing of such an aid is only to be indicated after consultation with the physician:

- Skin diseases and injuries to body parts to which the orthosis is to be fitted, especially in case of inflammatory symptoms. This applies to prominent scars featuring swelling, redness, and overheating.
- Lymph flow disorders – also unclear soft tissue swelling in the area of the upper extremities
- Sensitivity loss in the area of upper extremities

10 Instructions for Use

It is imperative that the orthosis be used in accordance with the indication. The acceptance and collaboration of the patient play an important role. It is necessary to explain to the patient how to use the orthosis properly!

It is essential to observe the instructions of the attending physician.

Improper alterations to the orthosis are not permitted.

It is recommended that you always wear a T-shirt or undershirt under the orthosis.

The daily duration of use and total period of application must be determined by the attending physician. If not prescribed otherwise by the physician, it is recommended to use the orthosis during the learning phase or while normal muscle strength is being gradually restored for 1 to 2 hours daily. In the course of the therapy, the time of use can be increased step by step based on the consultation with the physician. The orthosis must be re-adjusted by the orthotist if necessary.

If the patient is sitting or lying for a longer period, the orthosis should be removed. When the Back Plate (A) slips upwards while the orthosis is being worn, the patient should pull it downwards to its proper place.

The patient should be instructed to consult the physician immediately when he/she notices any unusual changes.

11 Care Instructions

When delivering the orthosis, instruct the patient also that only the textile components of the orthosis are washable (after removal of the Back Plate), not the Back Plate itself.

Wash the textile components at 30 °C (86 °F) with a mild detergent. The hook & loop closures and Plastic Buckle (G) must be closed while washing the orthosis. The velcroed Back Pad (B), Waist Strap (D) and Shoulder Straps (C) must be detached from the Back Plate (A). We recommend using a washing bag for washing the textiles. Do not expose to the sun. Do not use oven or radiator heat for drying. Air-dry. Do not dry clean or iron textile components.

The Back Plate (A) can be cleaned on the coated side with a damp cloth. After drying, velcro the Back Pad (B) to the Back Plate (A) again and attach the Shoulder Strap (C) and Waist Strap (D) to the Back Plate following the numbered sequence (1–8).

12 Restrictions of Use

The 50R20 Dorso Osteo Care is designed for use on **one** patient. This orthosis is designed for contact with intact skin only.

Attention!

Parts that must be shaped to the physical contours of a patient as well as parts that come into direct contact with the skin may represent functional or hygienic risks when used by another person.

Functional risks are for example:

- Formation of cracks in the material (may be invisible to the eye if the paint is undamaged)
- Material fatigue due to alternating shock loading
- Reduced elongation (strain hardening) after the components have been cold-formed several times.

The application period or wearing time of the orthosis depends on medical indication.

13 Liability

The manufacturer's warranty applies only if the product has been used under the conditions and for the purposes described. The manufacturer recommends that the product be used and maintained according to the instructions for use.

14 CE conformity

This product meets the requirements of the 93/42/EWG guidelines for medical products. This product has been classified as a class I product according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Otto Bock with sole responsibility according to appendix VII of the guidelines.

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Otto Bock HealthCare GmbH