
50R50 Lumbo Direxa, 50R51 Lumbo Direxa Women, 50R52 Lumbo Direxa High

DE Gebrauchsanweisung	3
EN Instructions for Use	8
FR Instructions d'utilisation	12
IT Istruzioni per l'uso	17
ES Instrucciones de uso	22
NL Gebruiksaanwijzing	27
PL Instrukcja użytkowania	32
HU Használati utasítás	37
CS Návod k použití	42
SK Návod na používanie	47
TR Kullanma talimatı	52
RU Руководство по применению	56
ZH 使用说明书	62

1


2

3

4

	Größe / Size	Taillenumfang* / Waist circumference	
		cm	inch
	XS	70 – 80	27.6 – 31.5
	S	80 – 90	31.5 – 35.4
	M	90 – 100	35.4 – 39.4
	L	100 – 110	39.4 – 43.3
	XL	110 – 120	43.3 – 47.2
	XXL	120 – 130	47.2 – 51.2

* ca. 2 cm oberhalb des Beckenkamms / approx. 2 cm / 0.8 inch above the iliac crest

Material	Gewebe / Fabric: PA, EL, Lycra®
	Korsettstäbe / Corset Splints: POM

1 Vorwort

INFORMATION

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INFORMATION

Date of the last update: 2012-05-24

- Please read this document carefully.
- Follow the safety instructions.

These instructions for use provide important information on the fitting and application of the 50R50 Lumbo Direxa, 50R51 Lumbo Direxa Women and 50R52 Lumbo Direxa High lumbar orthoses.

2 Application

2.1 Intended Use

The orthosis is intended **exclusively** for orthotic fittings of the trunk and **exclusively** for contact with intact skin.

The orthosis must be used according to the indications and following the instructions of a physician.

2.2 Indications

50R50 Lumbo Direxa, 50R51 Lumbo Direxa Women

- Moderate low back pain and sciatica
- Myostatic insufficiency/muscular imbalance of the lumbar trunk muscles
- Moderate degenerative changes of the lumbar spine (e.g. osteochondrosis, facet joint syndrome)
- Intervertebral disc protrusion L4-S1
- Stabilisation after discectomy L4-S1
- Moderate irritation of the sacroiliac joint

Indications must be determined by the physician.

50R52 Lumbo Direxa High

- Slight to moderate low back pain, sciatica, thoracic and lumbar ligament sprain
- Myostatic insufficiency/muscular imbalance of the lumbar and lower thoracic trunk muscles
- Moderate degenerative changes of the lumbar and lower thoracic spine (e.g. osteochondrosis, facet joint syndrome)
- Intervertebral disc protrusion L1-S1
- Stabilisation after discectomy L1-S1
- Moderate irritation of the sacroiliac joint

Indications must be determined by the physician.

2.3 Contraindications

2.3.1 Absolute Contraindications

Not known.

2.3.2 Relative contraindications

Consult the physician if any of the following occur: skin diseases/injuries, inflammation, raised scars with swelling, reddening, and warming in the area where the orthosis is fitted, tumours, lymphatic congestion – including unclear swelling of soft tissue distal to the device, and dysaesthesia in the legs.

2.4 Effects

The orthosis promotes sensorimotor functions and activates muscular stabilisation of the lumbar spine. It supports pain relief and relieves the lumbar spine.

3 Safety

3.1 Explanation of Warning Symbols

 CAUTION	Warning regarding possible risks of accident or injury
 NOTICE	Warning regarding possible technical damage

3.2 General Safety Instructions

CAUTION

Using the product without following the instructions for use

Deteriorating health condition and damage to the product because of failure to follow the safety information

- ▶ Observe the safety information in these instructions for use.
- ▶ Make the patient aware of all safety instructions marked with "**Inform the patient**".

CAUTION

Contact with heat, embers or fire

Injuries (e.g. burns) caused by melting of the material

- ▶ Keep the product away from open flame, embers and other sources of heat.
- ▶ **Inform the patient.**

CAUTION

Reuse on other patients and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ Only use the product for a single patient.
- ▶ Clean the product regularly.
- ▶ **Inform the patient.**

NOTICE

Contact with salves, lotions or other products that contain oils or acids

Insufficient stabilisation due to loss of material functionality

- ▶ Do not expose the product to salves, lotions or other products that contain oils or acids.
- ▶ **Inform the patient.**

4 Handling

INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

4.1 Selecting the size

- 1) Measure the waist circumference approx. 2 cm above the iliac crest.
- 2) Determine the orthosis size (see size table).

4.2 Fitting and Application

CAUTION

Direct skin contact with the product

Skin irritation due to chafing or sweat formation

- ▶ Do not wear this product directly on the skin.
- ▶ **Inform the patient.**

CAUTION

Incorrect or excessively tight application

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the orthosis is applied properly and fits correctly.
- ▶ **Inform the patient.**

NOTICE

Use of a worn or damaged product

Limited effectiveness

- ▶ Instruct the patient to inspect the product for proper functionality, wear and damage prior to each use.
- ▶ Advise the patient that the product must not be used if the product or part of the product exhibits signs of wear (e.g. cracks, deformation, poor fit) or damage.

A pad is also available in addition to the orthosis and can be ordered as a service part if necessary.

- 1) **Optional – install the pad:** Use the hook and loop closure to attach the pad to the inside of the orthosis so that it is at least 1 cm above the seat surface when seated.
- 2) **Applying** (see Fig. 1): Place the orthosis around the trunk. Place the fingers in the hand straps of the abdominal closure. Apply equal force to both sides to pull the orthosis forward.
- 3) **Hook and loop material** (see Fig. 2): Attach the right part to the top of the left part of the abdominal closure. Remove the fingers from the left and then from the right hand strap. Ensure that the hook and loop material sections lie flat against each other and are not in contact with clothing or skin.
- 4) **Increasing compression:** Pull the ends of the elastic straps forward (see Fig. 3) and attach them to the bottom closure (see Fig. 4).
- 5) **Final check:** The correct fit and proper positioning of the orthosis must be checked during handover to the patient. The orthosis is seated correctly when the lower edge of the orthosis extends just below the groin.
- 6) **Instruct the patient:** Explain to the patient how to don the orthosis and ensure that it is seated correctly. Practice donning and doffing of the orthosis with the patient.

4.3 Cleaning

Clean the orthosis regularly.

- 1) Close all hook and loop closures.
- 2) **Recommendation:** Use a laundry bag or net.
- 3) Wash the orthosis in warm water at 40 °C with standard mild detergent. Do not use fabric softener. Rinse thoroughly.
- 4) Allow to air dry. Do not expose to direct heat (e.g. from the sun, an oven or radiator).

5 Disposal

Dispose of the product in accordance with national regulations.

6 Legal Information

6.1 Liability

The manufacturer's warranty applies only if the device has been used under the conditions and for the purposes described. The manufacturer recommends that the device be used and maintained according to the instructions for use.

6.2 CE Conformity

This device meets the requirements of the 93 / 42 / EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to appendix VII of the guidelines.

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Otto Bock HealthCare GmbH