

50R59 Dorso Direxa Posture


DE Gebrauchsanweisung	4
EN Instructions for Use	8
FR Instructions d'utilisation	13
IT Istruzioni per l'uso	18
ES Instrucciones de uso	23
NL Gebruiksaanwijzing	28
PL Instrukcja użytkowania	33
HU Használati utasítás	38
CS Návod k použití	43
SK Návod na používanie	47
TR Kullanma talimatı	52
RU Руководство по применению	56
JA 取扱説明書	62
ZH 使用说明书	66

	Größe / Size	Taillenumfang* / Waist circumference	
		cm	inch
	XS	70 – 80	27.6 – 31.5
	S	80 – 90	31.5 – 35.4
	M	90 – 100	35.4 – 39.4
	L	100 – 110	39.4 – 43.3
	XL	110 – 120	43.3 – 47.2
	XXL	120 – 130	47.2 – 51.2

* ca. 2 cm oberhalb des Beckenkamms / approx. 2 cm/0.8 inch above the iliac crest

Material	Gewebe / Fabric : PA, EL, Lycra®
	Frotteebezüge / Terry Cloth Covers: Baumwolle / Cotton, PES
	Korsettstäbe / Corset Splints : Metall / Metal, POM
	Pelotte / Pad: TPE

1 Vorwort

INFORMATION

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INFORMATION

Date of the last update: 2012-05-24

► Please read this document carefully.

► Follow the safety instructions.

These instructions for use provide important information on the fitting and application of the 50R59 Dorso Direxa Posture thoracolumbar orthosis.

2 Application

2.1 Intended Use

The orthosis is intended **exclusively** for orthotic fittings of the trunk and **exclusively** for contact with intact skin.

The orthosis must be used according to the indications and following the instructions of a physician.

2.2 Indications

- Poor body posture (hunchback without fused vertebrae, Scheuermann's disease)
- Osteoporosis of the thoracic and lumbar spine
- Aftercare of osteoporotic fractures of the thoracic and/or lumbar spine
- Back pain (lumbar and lower thoracic spine), sciatica, muscular imbalance of the thoracic and lumbar spine, thoracic and lumbar spine ligament sprain
- Slight degenerative changes of the thoracic and lumbar spine (e.g. osteochondrosis, facet joint syndrome)

Indications must be determined by the physician.

2.3 Contraindications

2.3.1 Absolute Contraindications

Not known.

2.3.2 Relative contraindications

The following indications require consultation with a physician: skin diseases/injuries, inflammation, raised scars with swelling, reddening and warming of the treated body areas, tumours, lymphatic flow disorders - including unclear soft tissue swelling in areas distant to where the device is applied, sensitivity disorders in the arms and legs.

2.4 Effects

The orthosis supports muscular stabilisation of the thoracic spine and straightens it actively. It promotes sensorimotor functions, supports pain relief and relieves the thoracic spine through active straightening and muscular stabilisation and the lumbar spine by increasing intra-abdominal pressure and providing external stabilisation.

3 Safety

3.1 Explanation of Warning Symbols

⚠ CAUTION Warning regarding possible risks of accident or injury

NOTICE Warning regarding possible technical damage

3.2 General Safety Instructions

⚠ CAUTION

Using the product without following the instructions for use

Deteriorating health condition and damage to the product because of failure to follow the safety information

- ▶ Observe the safety information in these instructions for use.
- ▶ Make the patient aware of all safety instructions marked with **"Inform the patient"**.

⚠ CAUTION

Contact with heat, embers or fire

Injuries (e.g. burns) caused by melting of the material

- ▶ Keep the product away from open flame, embers and other sources of heat.
- ▶ **Inform the patient.**

⚠ CAUTION

Reuse on other patients and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ Only use the product for a single patient.
- ▶ Clean the product regularly.
- ▶ **Inform the patient.**

NOTICE

Contact with salves, lotions or other products that contain oils or acids

Insufficient stabilisation due to loss of material functionality

- ▶ Do not expose the product to salves, lotions or other products that contain oils or acids.
- ▶ **Inform the patient.**

4 Handling

INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

4.1 Selecting the size

- 1) Measure the waist circumference approx. 2 cm above the iliac crest.
- 2) Determine the orthosis size (see size table).

4.2 Fitting and Application

CAUTION

Direct skin contact with the product

Skin irritation due to chafing or sweat formation

- ▶ Do not wear this product directly on the skin.
- ▶ **Inform the patient.**

CAUTION

Incorrect or excessively tight application

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the orthosis is applied properly and fits correctly.
- ▶ **Inform the patient.**

NOTICE

Use of a worn or damaged product

Limited effectiveness

- ▶ Instruct the patient to inspect the product for proper functionality, wear and damage prior to each use.
- ▶ Advise the patient that the product must not be used if the product or part of the product exhibits signs of wear (e.g. cracks, deformation, poor fit) or damage.

- 1) **Fitting the metal rods** (see Fig. 1): Remove the 2 metal rods from the pockets and fit them to the spine individually. Put the metal rods back in the pockets.
- 2) **Optional – install the pad** (see Fig. 2): Use the hook and loop closure to attach the pad to the inside of the orthosis so that it is at least 1 cm above the seat surface when seated.
- 3) **Application:** Place the upper part of orthosis over the shoulders like a backpack (see Fig. 3). Place the lower part of the orthosis around the trunk. Place the fingers in the hand straps of the abdominal closure. Apply equal force to both sides to pull the orthosis forward (see Fig. 4).
- 4) **Hook and loop material:** Attach the right part to the top of the left part of the abdominal closure (see Fig. 5). Remove the fingers from the left and then from the right hand strap. Ensure that the hook and loop material sections lie flat against each other and are not in contact with clothing or skin. Route the ends of the shoulder strap to the front around the chest (see Fig. 6) and secure the hook and loop closure below the chest (see Fig. 7).
- 5) **Increasing compression:** Pull the ends of the elastic straps forward (see Fig. 8) and attach them to the bottom closure (see Fig. 9).
- 6) **Final check:** The correct fit and proper positioning of the orthosis must be checked during handover to the patient. The orthosis is seated correctly when the lower edge of the orthosis extends just below the groin.
- 7) **Instruct the patient:** Explain to the patient how to don the orthosis and ensure that it is seated correctly. Practice donning and doffing of the orthosis with the patient.

4.3 Cleaning

Clean the orthosis regularly.

- 1) Remove the metal rods. Close all hook and loop closures.
- 2) **Recommendation:** Use a laundry bag or net.
- 3) Wash the orthosis in warm water at 40 °C with standard mild detergent. Do not use fabric softener. Rinse thoroughly.
- 4) Allow to air dry. Do not expose to direct heat (e.g. from the sun, an oven or radiator).

5 Disposal

Dispose of the product in accordance with national regulations.

6 Legal Information

6.1 Liability

The manufacturer's warranty applies only if the device has been used under the conditions and for the purposes described. The manufacturer recommends that the device be used and maintained according to the instructions for use.

6.2 CE Conformity

This device meets the requirements of the 93 / 42 / EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to appendix VII of the guidelines.

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Otto Bock HealthCare GmbH