

50R49 Dorso Carezza Posture

DE Gebrauchsanweisung	3
EN Instructions for Use	7
FR Instructions d'utilisation	11
IT Istruzioni per l'uso	15
ES Instrucciones de uso	19
NL Gebruiksaanwijzing	23
PL Instrukcja użytkowania	28
HU Használati utasítás	32
CS Návod k použití	36
SK Návod na používanie	40
TR Kullanma talimatı	44
RU Руководство по применению	48
JA 取扱説明書	52
ZH 使用说明书	56

	Größe / Size	Taillenumfang* / Waist circumference	
		cm	inch
	XXS	60 – 70	23.6 – 27.6
	XS	70 – 80	27.6 – 31.5
	S	80 – 90	31.5 – 35.4
	M	90 – 100	35.4 – 39.4

* ca. 2 cm oberhalb des Beckenkamms / approx. 2 cm/0.8 inch above the iliac crest

Material	Gewebe / Fabric: PA, EL, Lycra®
	Korsettstäbe / Corset Splints: POM

1 Vorwort

INFORMATION

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INFORMATION

Date of the last update: 2012-05-24

- ▶ Please read this document carefully.
- ▶ Follow the safety instructions.

These instructions for use provide important information on the fitting and application of the 50R49 Dorso Carezza Posture thoracic spine orthosis.

2 Application

2.1 Intended Use

The orthosis is intended **exclusively** for orthotic fittings of the trunk and **exclusively** for contact with intact skin.

The orthosis must be used according to the indications and following the instructions of a physician.

2.2 Indications

- Poor body posture (hunchback without fused vertebrae, Scheuermann's disease)
- Back pain (thoracic spine)

Indications must be determined by the physician.

2.3 Contraindications

2.3.1 Absolute Contraindications

Not known.

2.3.2 Relative contraindications

Consult the physician if any of the following occur: skin diseases/injuries, inflammation, raised scars with swelling, reddening, and warming in the area where the orthosis is fitted, tumours, lymphatic congestion – including unclear swelling of soft tissue distal to the device, and dysaesthesia in the legs.

2.4 Effects

The orthosis supports muscular stabilisation of the thoracic spine and straightens it actively.

3 Safety

3.1 Explanation of Warning Symbols

CAUTION Warning regarding possible risks of accident or injury

NOTICE Warning regarding possible technical damage

3.2 General Safety Instructions

Using the product without following the instructions for use

Deteriorating health condition and damage to the product because of failure to follow the safety information

- ▶ Observe the safety information in these instructions for use.
- ▶ Make the patient aware of all safety instructions marked with **"Inform the patient"**.

Contact with heat, embers or fire

Injuries (e.g. burns) caused by melting of the material

- ▶ Keep the product away from open flame, embers and other sources of heat.
- ▶ **Inform the patient.**

Reuse on other patients and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ Only use the product for a single patient.
- ▶ Clean the product regularly.
- ▶ **Inform the patient.**

NOTICE

Contact with salves, lotions or other products that contain oils or acids

Insufficient stabilisation due to loss of material functionality

- ▶ Do not expose the product to salves, lotions or other products that contain oils or acids.
- ▶ **Inform the patient.**

4 Handling

INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

4.1 Selecting the size

- 1) Measure the waist circumference approx. 2 cm above the iliac crest.
- 2) Determine the orthosis size (see size table).

4.2 Fitting and Application

⚠ CAUTION

Direct skin contact with the product

Skin irritation due to chafing or sweat formation

- ▶ Do not wear this product directly on the skin.
- ▶ **Inform the patient.**

⚠ CAUTION

Incorrect or excessively tight application

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the orthosis is applied properly and fits correctly.
- ▶ **Inform the patient.**

NOTICE

Use of a worn or damaged product

Limited effectiveness

- ▶ Instruct the patient to inspect the product for proper functionality, wear and damage prior to each use.
- ▶ Advise the patient that the product must not be used if the product or part of the product exhibits signs of wear (e.g. cracks, deformation, poor fit) or damage.

- 1) **Application** (see Fig. 1): Place the orthosis over the shoulders like a backpack.
- 2) **Hook and loop material:** Route the ends of the lower strap to the front (see Fig. 2) and secure the hook and loop closure (see Fig. 3). Route the ends of the shoulder strap to the front around the chest and secure the hook and loop closure below the chest (see Fig. 4).
- 3) **Final check:** The correct fit and proper positioning of the orthosis must be checked during handover to the patient.
- 4) **Instruct the patient:** Explain to the patient how to don the orthosis and ensure that it is seated correctly. Practice donning and doffing of the orthosis with the patient.

4.3 Cleaning

Clean the orthosis regularly.

- 1) Close all hook and loop closures.
- 2) **Recommendation:** Use a laundry bag or net.
- 3) Wash the orthosis in warm water at 40 °C with standard mild detergent. Do not use fabric softener. Rinse thoroughly.
- 4) Allow to air dry. Do not expose to direct heat (e.g. from the sun, an oven or radiator).

5 Disposal

Dispose of the product in accordance with national regulations.

6 Legal Information

6.1 Liability

The manufacturer's warranty applies only if the device has been used under the conditions and for the purposes described. The manufacturer recommends that the device be used and maintained according to the instructions for use.

6.2 CE Conformity

This device meets the requirements of the 93 / 42 / EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to appendix VII of the guidelines.

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Otto Bock HealthCare GmbH